
User manual

Iriscope

Version 1 – 2023-04

Languages

Pages

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Contents

1	Important Information.....	32
1.1	Intended use.....	32
1.2	User qualifications.....	32
1.3	Indications for use.....	32
1.4	Precautions for use.....	32
2	Description.....	33
3	Before use.....	33
3.1	Storage.....	33
3.2	Compatibility with endoscopes/endoscopic tools.....	33
4	Use.....	34
4.1	Preloading the catheter.....	34
4.2	Inserting the catheter in the endoscope.....	35
4.3	Handling the catheter.....	35
4.4	Use with endoscopic tools.....	36
4.5	Removal of the catheter.....	36
5	Troubleshooting.....	36
6	Disposal of the catheter.....	37
7	Technical specifications.....	37
7.1	Specifications.....	37
7.2	Applicable standards.....	37
8	Symbols used.....	37

1 Important Information

These operating instructions must be read **carefully** prior to using the *Iriscope GS* catheter. They contain the information required to use it safely.

This user manual describes the *Iriscope GS* catheter, its operation and its safety precautions. It is not intended to describe clinical procedures. It is thus imperative that users have been trained in endoscopic techniques to use the *Iriscope GS*, and be aware of all precautions, warnings, indications et contra-indications set forth in this user manual.

All serious incidents related to the *Iriscope GS* catheter must be reported immediately to Lys Medical (vigilance@lysmedical.com) and to the responsible authority in the country where you are based.

Failure to follow the operating instructions could lead to injury to the patient and/or damage to the device. Lys Medical shall not be held liable in the event of improper use.

1.1 Intended use

The *Iriscope GS* catheter is intended for use through interventional bronchoscopes to extend the operating channel and thus serve as a support for compatible endoscopic tools (e.g. *Iriscope* probe, biopsy forceps, etc.). It is compatible with Olympus (190 series) interventional bronchoscopes, or equivalent, with a 2.8 mm operating channel diameter or larger.

1.2 User qualifications

This device is intended for use by a physician specializing in bronchial endoscopy.

These operating instructions only deal with the operation of the *Iriscope GS* catheter and its use, and endoscopic treatment and techniques are not the subject of this user manual.

1.3 Indications for use

Iriscope GS is a single-use navigation catheter. It must be discarded after use.

Endoscopic tools (brush, biopsy forceps, probes) can be inserted into the navigation catheter. Tools up to a diameter of 2.0 mm can be inserted. However, it is recommended to test the insertion of tools in the guide beforehand (see section 3.2).

The catheter is reserved for use in adults, for the exploration of the respiratory tract, as far as the sub-segmental bronchi.

1.4 Precautions for use

It is important to heed these warnings and precautions to ensure patient safety and the service life of the device.

Also take care to carefully read the precautions regarding catheter storage in section 3.1.

Warnings and precautions

- The catheter must only be used by physicians trained in clinical endoscopic procedures and techniques.
- Before use, check that the packaging is not damaged. Never use the catheter if the packaging is damaged.
- Always check the compatibility of the *Iriscope GS* catheter with the endoscope operating channel in advance, as well as the compatibility of endoscopic tools with *Iriscope GS*.
- Always preload the catheter with an endoscopic tool before use. Inserting the catheter empty could damage the catheter, endoscope, or the patient's anatomical pathways.
- Never force the insertion or removal of the *Iriscope GS* catheter in the endoscope operating channel.
- Never force the insertion or removal of endoscopic tools in the *Iriscope GS* catheter.
- In case of malfunction during the procedure, halt the procedure and remove the *Iriscope GS* catheter from the patient.
- Be careful not to bend the catheter and not to drop the handle as this could permanently damage it.
- Do not use the device if the expiry date has passed. After this date, the cleanliness of the product is no longer guaranteed.

2 Description

Iriscope GS is a navigation catheter supporting *Iriscope* probes or any other endoscopic tool.


Figure 1 : *Iriscope GS* catheter and handle, with (1) flush access, (2) rotation handle controlling the orientation of the catheter, (3) attachment of the catheter to the endoscope, (4) attachment to the operator channel of the endoscope, and (5) navigation catheter.

This catheter is intended for exploring the airways for a maximum duration of 60 minutes.

3 Before use

3.1 Storage

⚠ Warnings and precautions

- Be sure to comply with the storage conditions set out in this manual as improper storage could lead to an infection control risk and/or damage to the device.
- Ensure that the catheter is kept upright and flat during transport and storage. Non-compliant storage could permanently damage its shape by inducing catheter deformation over time and reduced transmission of rotation.
- Take care not to leave the catheter exposed to direct sunlight (UV) as this could speed up its deterioration and greatly reduce the lifetime of the catheter, or even render it unusable.

The catheter in its original packaging must be kept in a dry place away from moisture and light in a clean room, at standard temperature and pressure.

3.2 Compatibility with endoscopes/endoscopic tools

It is important to check the compatibility between the *Iriscope GS* catheter and the operating channel of the endoscope and endoscopic tools used:

- The outer diameter of the catheter is 2.7 mm, compatible with a 2.8 mm endoscope operating channel.
- The inner diameter of the catheter is 2.1 mm, compatible with endoscopic tools with an outer diameter of up to 2.0 mm.

However, it cannot be guaranteed that instruments selected on the basis of the minimum or maximum recommended diameter will be compatible with the catheter. It is therefore advisable to perform a preliminary insertion test with the tools you wish to use to ensure their mechanical compatibility.

To do so, insert the catheter into the endoscope and check that insertion takes place without forcing, and that the catheter slides in correctly without catching. The use of medical lubricant can help with catheter insertion. Likewise, insert the endoscopic tools into the catheter and check that insertion takes place without forcing, and that the tools slide correctly without catching. The use of medical lubricant can help with the insertion of tools.


Incompatibility could result in damage to the *Iriscope GS* catheter and/or injury to the patient.

4 Use


4.1 Preloading the catheter

Before inserting the catheter into the endoscope:

1. Insert the compatible endoscopic tool (*Iriscope* probe, forceps or other) into the *Iriscope* GS catheter channel. Ensure that the tip of the tool is flush with the distal end of the *Iriscope* GS catheter.

2. Tighten the screw cap of the Tuohy-Borst adapter for flushing to hold the endoscopic tool in position. Flushing with a syringe allows the seal to be tested.

3. Check that the sliding lock is in the open position (Figure 2Figure 2).



Figure 2 : Open position (left) and closed position (right) of the endoscope handle locking mechanism.

4.2 Inserting the catheter in the endoscope

1. Check that the operating channel tip is freely accessible (no sealing valve).

2. Insert the catheter tip into the endoscope operating channel. Be careful when pushing the catheter to avoid bending it, which would damage it irreversibly. Rotating while inserting the catheter makes it easier to insert it into the endoscope operating channel. The use of medical lubricant also facilitates its insertion. If resistance is encountered, do not force it – remove the catheter to check where the resistance is coming from.

Once the catheter is inserted into the endoscope operating channel, the handle can be attached to the operating channel. Start by sliding the handle downwards on to the operating channel pin (Figure 3Figure 3-1), inserting it to the end (Figure 3Figure 3-2). Move the sliding lock to the closed position (Figure 3Figure 3-3). Check that the handle is securely attached to the endoscope (no vertical movements or rotation).


Figure 3 : Inserting the handle on to the endoscope operating channel and locking it.

4.3 Handling the catheter

The catheter does not have an angulation system but is preformed.

The handle facilitates the movement of the navigation catheter. The catheter is secured to the handle. It is possible to rotate the catheter by the handle, as well as perform translational movement over a stroke of 18 cm (Figure 4Figure 4). Be careful not to force if resistance is encountered, as this could irreversibly damage the equipment.


Figure 4 : Handling the Iriscope GS catheter with the rotation handle and with the telescopic handle over an 18 cm stroke.

4.4 Use with endoscopic tools

Once in position or during the procedure, the tools pre-loaded in the *Iriscope GS* catheter can be removed. Different tools can be inserted or removed during the procedure as long as they are compatible with the dimensions of the guide catheter.

In case of pronounced bending of the navigation catheter, a risk of kinking may arise. Bending the catheter to a radius of curvature of less than 10 mm is not recommended. In case of kinking, it may not be possible to insert other tools. Likewise, in the event of kinking, it may be difficult to remove tools. It is recommended to remove the guide catheter from the endoscope to de-kink the catheter in order to remove trapped tools. Dispose of the kinked *Iriscope GS* catheter.

Take care not to force when inserting or removing endoscopic tools if resistance is encountered, as this could irreversibly damage the equipment.

The Tuohy-Borst adapter (Figure 5) allows water or saline solution to be flushed around the tool inserted into the catheter during use. Once in position, tighten the clamping ring (without forcing to avoid damaging the equipment). Use a syringe with a tip compatible with "Luer" type connectors. Check that there is no proximal flow.



Figure 5 : Tuohy-Borst adapter with (1) the clamping ring that allows the tightness around the probe and the insertion friction to be adjusted and (2) access for the syringe.

4.5 Removal of the catheter

Place the sliding lock in the open position. Pull the handle upwards while keeping the endoscope in position. Carefully remove the catheter from the endoscope.

In the event of blockage, rotate or oscillate to try to release the catheter.

5 Troubleshooting

Faults	Possible causes	Solutions
The handle does not attach properly to the endoscope.	The sliding lock is in the closed position.	Place the sliding lock in the open position (check the lower part of the handle: a circular opening must be visible).
	A biopsy valve is placed on the operating channel.	Remove the valve to free the operating channel.

Faults	Possible causes	Solutions
	The endoscope is not compatible.	Use the catheter alone or without attaching the handle. Alternatively, use a compatible endoscope (Olympus 190 type).
The passage of an endoscopic tool is blocked.	The catheter is kinked.	Replace the catheter.
	The diameter of the tool is not compatible.	Use another tool with a compatible diameter.

6 Disposal of the catheter

This is a single-use catheter. After use, dispose of it in accordance with your institution's guidelines and all applicable regulations.

7 Technical specifications

7.1 Specifications

Iriscope GS catheter	<i>Iriscope GS 90</i>	<i>Iriscope GS 180</i>
Maximum outer diameter of the inserted part [mm]	2.7	2.7
Nominal inner diameter [mm]	2.05	2.05
Angle of curvature (preformed part) [°]	90	180
Working length [cm]	80	80
Overall length [cm]	105	105
Minimum radius of curvature [mm]	10	10
Compatibility		
Maximum diameter of endoscopic tools [mm]	2.0	2.0
Minimum diameter of endoscope operating channel [mm]	2.8	2.8
Storage and transport		<i>Iriscope GS</i>
Storage temperature [°C]	0 – 40	
Transport temperature [°C]	-20 – 60	
Storage relative humidity [%]	30 – 80	
Transport relative humidity [%]	10 – 90	
Atmospheric pressure [kPa]	80 – 110	

7.2 Applicable standards

The *Iriscope GS* catheter operates in accordance with European regulation *MDR 2017/745* and with the following standards:

- ✓ ISO 10993-1: Biological evaluation of medical devices — Part 1: Evaluation and testing within a risk management process

8 Symbols used

Symbols	Meaning
	Medical device.
	Unique device identification.
	Catalogue reference.
	CE marking. The device complies with European regulation <i>MDR 2017/745</i> .
	Medical device manufacturer.
	Batch number.

Symbols	Meaning
	Use-by date.
	Non-sterile device.
	Do not reuse.
	Keep out of direct sunlight.
	See instructions for use.
	Do not use if the packaging is damaged.
	Protect from moisture.
	Fragile, handle with care.
	Catheter preforming angle.
	<i>Iriscope GS</i> inner diameter.
	<i>Iriscope GS</i> outer diameter.
	Working length, or overall length.

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