User manual

# Iriscope<sup>®</sup>

Version 5 – 2023-06

Languages	Pages
Français (Instructions d'utilisation).....	05-17
Nederlands (Gebruiksaanwijzing).....	18-30
Deutsch (Bedienungsanleitung).....	31-43
English (Operating instructions).....	44-55
Español (Instrucciones de uso).....	56-68
Italiano (Istruzioni d'uso).....	69-81

## Contenu

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## Contents

1.1	Intended use .....	45
1.2	User qualifications .....	45
1.3	Indications for use .....	45
1.4	Precautions for use .....	45
2	Description .....	46
3	New probe and precautions before first use .....	47
3.1	Cleaning before first use .....	47
3.2	Storage .....	47
3.3	Compatibility with endoscopes/guide sleeves .....	47
4	Use .....	47
4.1	Preparation and inspection .....	47
4.1.1	Visual inspection of the probe .....	47
4.1.2	Connecting the probe and inspecting the image .....	48
4.2	Operation .....	48
4.2.1	Probe insertion .....	48
4.2.2	Probe handling .....	48
4.2.3	Probe removal .....	48
4.3	After use .....	48
4.4	Reprocessing of the probe .....	48
4.4.1	Manual Reprocessing – Quaternary ammoniums .....	49
4.4.2	Automatic Reprocessing .....	49
4.4.3	Hydrogen peroxide sterilisation .....	50
4.5	Cleaning the connector .....	50
4.5.1	Cleaning electrical contacts .....	50
4.5.2	Cleaning the ferrule (white cylinder) .....	50
4.6	Storing the probe .....	51
4.7	Troubleshooting .....	51
5	Disposal of the probe .....	51
6	Technical specifications .....	51
6.1	Specifications .....	51
6.2	Applicable standards .....	52
7	Symbols used .....	52
8	Appendix .....	54
A.	Instructions for manual reprocessing .....	54
B.	Instructions for automatic reprocessing - .....	55

# 1 Important Information

These operating instructions must be read **carefully** prior to using the *Iriscope* probe. They contain the information required to use it safely.

This user manual describes the *Iriscope* probe, its operation and its safety precautions. It is not intended to describe clinical procedures. It is thus imperative that users should be trained in endoscopic techniques to use the *Iriscope*, and be aware of all precautions, warnings, indications and contra-indications set forth in this user manual.

All serious incidents related to the *Iriscope* probe must be reported immediately to Lys Medical ([vigilance@lysmedical.com](mailto:vigilance@lysmedical.com)) and to the responsible authority in the country where you are established.

Failure to follow the operating instructions could lead to injury to the patient and/or damage to the device. Lys Medical shall not be held liable in the event of improper use.

## 1.1 Intended use

The *Iriscope* probe (including *Iriscope17*, *Iriscope18*, and *Iriscope20*) is designed to visualise the respiratory tract and the upper gastrointestinal tract. This probe can only be used in combination with the *Iristar* control box. It is intended to be used during an endoscopic procedure, with an endoscope or an endoscopic accessory.

*Iriscope* is a very small video endoscope enabling visualising in small tracts.

## 1.2 User qualifications

The device may only be used in a hospital by, or under the supervision of a specialised endoscopy physician. The user's qualifications to perform an endoscopy must follow the official guidelines, standards or laws established by the local or international authorities, or by a medical association.

The user is therefore qualified to perform an endoscopy including the risks or complications that may arise during endoscopy. These operating instructions only deal with the operation of the *Iriscope* probe and its use, and endoscopy treatment techniques are not the purpose of this user manual.

## 1.3 Indications for use

*Iriscope* is a reusable endoscopic probe but has a limited lifetime. A visual inspection of the probe must be carried out before each use to ensure its good condition. Special attention should be paid to the catheter. See section 4.1, and more particularly 4.1.1, for inspecting the probe. In case of damage the probe cannot be reused and must be duly disposed of.

Due to its reuse, the probe must be cleaned and disinfected after each use following the method presented in this user manual (see section 4.4).

The *Iriscope17* probe can be inserted into any guide sleeve with an operating channel whose diameter is greater than or equal to 1.9 mm. The *Iriscope18* probe can be inserted into any guide sleeve with an operating channel whose diameter is greater than or equal to 2.0 mm. The *Iriscope20* probe can be inserted into any guide sleeve with an operating channel whose diameter is greater than or equal to 2.1 mm. However, it is recommended to insert the *Iriscope* probe into the guide sleeve before use on the patient (see section 3.3).

The probe is reserved for use in adults, for the exploration of the respiratory tract, as far as the subsegmental bronchi, and of the gastrointestinal tract, as far as the upper duodenum.

## 1.4 Precautions for use

It is important to heed these warnings and precautions to ensure patient safety and the service life of the device.

Care should also be taken to carefully read the precautions regarding probe storage and reprocessing, which can be found directly in the sections concerning these points (see section 3.2 and section 4.4, respectively).

### Warnings and precautions

- Before the first use, check that the packaging is not damaged. Do not use the probe if the packaging is damaged.
- Before and after use, perform a visual check that the *Iriscope* probe is not damaged (see section 4.1). Never use the probe if it is damaged.
- Before the first use, clean and disinfect the probe following the instructions in this manual.
- The probe can only be used by physicians trained in clinical endoscopy procedures and techniques.
- It cannot be guaranteed that the instruments selected on the basis of the minimum recommended diameter will be compatible with the probe. Always check the compatibility between the *Iriscope* probe and the operating channel of the endoscope or guide sleeve before use in the patient.

- Never force the insertion or removal of the *Iriscope* probe in the operating channel.
- Do not use the *Iriscope* probe if the connector plastic protection is broken.
- Do not use the *Iriscope* probe in an MRI environment.
- Do not use the *Iriscope* probe during defibrillation.
- Do not use the *Iriscope* probe with high-frequency surgical equipment or with equipment using laser or ultrasound.
- Do not use the *Iriscope* probe while administering a highly inflammable anaesthetic gas to the patient. This could cause injury to the patient.
- High-intensity light is emitted from the tip of the probe, which can generate high temperatures at its front (up to 50°C). To reduce the risk of burns, avoid putting the tip of the probe in prolonged contact with the same area of the mucosa.
- Do not apply pressure with the distal end of the probe to prevent heat transfer to the mucosa.
- Never target the eyes when the probe is connected and lit as it may dazzle.
- In case of malfunction during the procedure, halt the procedure and remove the *Iriscope* probe from the patient.
- Never disconnect the probe from its control box while it is in the patient.
- Never clean the connector when the probe is in the patient.
- Never hold the probe by its catheter letting the connector hang from its end during transport or handling as this could permanently damage it.
- Be careful not to bend the probe and not to drop the connector as this could permanently damage it.
- Improper handling, such as touching the reprocessed probe with contaminated gloves, placing the reprocessed device on a contaminated support or surface, allowing the device to touch the ground, etc., will recontaminate the device.

## 2 Description

*Iriscope* is a very small video endoscope with a lighting system. This probe must be connected to the *Iristar* control box to acquire images and emit light. The *Iriscope* probe is an applied part of the *Iristar* control box. Please refer to the operating instructions of this control box for further information.


Figure 1 - *Iriscope* probe

This probe is intended to explore the natural tracts, including the respiratory tract and the upper gastrointestinal tract, for a maximum time of 60 minutes. It has been designed for use by endoscopic practitioners (ENT, gastroenterologists, pneumologists) or any physician trained to perform endoscopies. It can be used in any patient who may require an endoscopy.

The probe has a protective cap for its connector. It must be placed on the control box connector as soon as it is not in use in order to prevent it from getting dirty.

## 3 New probe and precautions before first use

### 3.1 Cleaning before first use

The new probe must be cleaned and disinfected before its first use. Please see section 4.4 to perform cleaning/disinfection.

### 3.2 Storage

#### Warnings and precautions – Stockage

- Be sure to comply with the storage conditions set out in this manual as improper storage could lead to an infection control risk and/or device degradation.
- Make sure that the probe winding is loose, with a curvature radius greater than 10 cm, during transport and storage. Coiling with too small a curvature radius could cause permanent damage to its shape leading the catheter to deform over time.
- Take care not to leave the catheter exposed to direct sunlight (UV) as this could speed up its deterioration and greatly reduce the lifetime of the probe, or even render it unusable.

The new probe in its original packaging must be kept in a dry place away from moisture and light in a clean room, at standard temperature and pressure.

After opening the packaging and between each use, the probe must be kept in a clean, well-ventilated room, protected from moisture and light, and at standard temperature and pressure.

The probe can be stored by being loosely rolled (curvature radius approximately 10 cm) and placed in a clean tray.

If it is possible to store the probes lengthwise thus bending them less, then this is preferable.

### 3.3 Compatibility with endoscopes/guide sleeves

It is important to check the compatibility between the *Iriscope* probe and the operating channel of the endoscope or guide sleeve used:

- The *Iriscope17* probe is intended for insertion into a guide sleeve with a minimum operating channel of 1.9 mm,
- *Iriscope18* with a minimum operating channel of 2.0 mm, and
- *Iriscope20* with a minimum operating channel of 2.1 mm.

However, it cannot be guaranteed that the instruments selected on the basis of the minimum recommended diameter will be compatible with the probe. It is therefore advisable to perform a preliminary insertion test with the guide sleeve you wish to use to ensure their mechanical compatibility.

To do so, insert the probe into the guide sleeve and check that insertion takes place without forcing, and that the probe slides in correctly without blocking.

Incompatibility could result in damage to the *Iriscope* probe and/or injury to the patient.

## 4 Use

This section only presents the use of the *Iriscope* probe, the instructions for using the *Iristar* control box are presented in a user manual of its own.

### 4.1 Preparation and inspection

#### 4.1.1 Visual inspection of the probe

1. Visually inspect the condition of the probe. Two things must be checked:
  - Check the catheter for damage, signs of deterioration, protrusions.
  - It is important that the metal filaments present in the catheter do not protrude. Slide your hand along the catheter to check this point; the catheter should be uniform.
  - Check that the probe head is secured to the catheter and the connector.
  - Check that the connector plastic protection is intact. It must be free of any cracks.

In case of damage, the probe cannot be reused and must be disposed of.

2. Remove the connector cap. Make sure that the inside of the connector is free of residual dirt and that it is dry. If the ferrule (white cylinder) looks dirty, clean it (see section 4.5.2).

#### 4.1.2 Connecting the probe and inspecting the image

1. Connect the probe to the control box by inserting the probe connector into the ad hoc control box connector by aligning them.
2. Check that an image appears on the screen connected to the control box. Point the probe tip at an object to check the quality of the image (e.g. your hand).
3. If the quality of the image is not good enough, check that the probe tip is clean and clean it with a sterile cloth if necessary.
4. Check the brightness; if it is not sufficient or lower than usual, check the cleanliness of the *Iristar* control box connectors and of the *Iriscope* probe and clean them if necessary (see section 4.5).

## 4.2 Operation

### 4.2.1 Probe insertion

Insert the probe tip into the endoscope or guide sleeve used with it. Be careful when pushing the probe to avoid bending its catheter, which would cause irreversible damage. Do not push too hard if there is any resistance.

If the image becomes cloudy, clean the tip after removing it.

### 4.2.2 Probe handling

The probe does not have a bending capability, so navigation is carried out with the guide sleeve or the endoscope into which it is inserted.

The probe has a reinforced catheter to allow it to rotate. To do so, gently pinch the probe between the fingers and rotate it. Be careful not to force it if it offers too much resistance; it could cause irreversible damage.

### 4.2.3 Probe removal

Carefully remove the probe from the guide sleeve or endoscope. Since the head of the probe is slightly wider than the catheter, be careful when retracting the catheter into the operating channel as it may get slightly stuck.

In the event of a blockage, perform a back-and-forth and/or swinging movement to try to release it. If this blockage persists, remove the guide sleeve or endoscope from the patient before removing the probe.

Do not remove it too quickly in the event of resistance and do not force it. This could distort the probe lengthwise and damage the lighting system.

## 4.3 After use

Detach the connector from the probe by pulling on the embossed part. Take care not to dirty it and immediately replace the watertight cap. Also put the watertight cap into the control box connector.

Clean the probe immediately with a cleaning product wipe (see section 4.3). Place it rolled up in a suitable tray to protect it and transport it to the cleaning area.

Then clean it in a suitable area according to the protocols in force in the place of use (see section 4.4).

## 4.4 Reprocessing of the probe

Different reprocessing methods are proposed in this manual. These have been validated for reprocessing the probe.

You can use other compatible methods, but please make sure that you use products certified for cleaning and disinfecting light-free thermosensitive endoscopic devices and follow the manufacturer's instructions. If you have any questions regarding the compatibility of a cleaning agent with the probe, please contact us.

The proposed methods are described in the following sections. There are two disinfection methods and one sterilisation method. Please note that the sterilisation method has a strong impact on the lifetime of the product.

### Warnings and precautions for reprocessing

- The probe is reusable so ensure it is cleaned and disinfected after each use according to the instructions herein.
- Be sure to follow the reprocessing instructions specified by the manufacturer, as well as those in this manual. Following the manufacturer's instructions ensures that the device is clean. Failure to do so could result in damage to the device
- Make sure that the probe connector is fitted with its cap before cleaning, otherwise it could damage the probe.
- The probe should be cleaned after each use and before each disinfection. If the probe is not properly cleaned, disinfection will not be effective.

- The probe must be pre-cleaned immediately after each procedure, in the endoscopy suite. Otherwise, residual organic debris will begin to solidify and this will hinder the effectiveness of cleaning and disinfection.
- An insufficiently cleaned/disinfected probe may pose an infection control risk to patients and/or operators who come into contact with it.
- If the disinfectant solution is reused, check its effectiveness with a test strip as recommended by the disinfectant manufacturer before use.
- Do not underdose the amount of detergent as reprocessing the probe may not be effective.
- Never reuse rinse water.
- Never exceed a temperature of 60°C as this could cause irreversible damage to the *Iriscope* probe as it is heat sensitive.
- Before using the probe, check that it has been properly disinfected and stored. If there is the slightest doubt, clean and disinfect the probe before use in a patient.
- Always remove contaminated personal protective equipment before leaving the reprocessing area to prevent contamination from spreading.
- Prions, the pathogen of Creutzfeldt-Jakob disease (CJD), cannot be destroyed or inactivated by cleaning, disinfection and sterilisation methods. When the probe is used on patients with CJD or a variant of Creutzfeldt-Jakob Disease (vCJD), be sure to use it only for these patients or dispose of it appropriately immediately after use to prevent the use of exposed apparatus on other patients. For methods for treating CJD; follow the respective guidelines of your country.

#### 4.4.1 Manual Reprocessing – Quaternary ammoniums

We propose manual reprocessing by soaking the probe in a quaternary ammonium solution, which is divided into three steps:

- Pre-cleaning: carried out directly after the procedure using a wipe (e.g. Umonium Neutralis Tissues from Huckert's International) to remove larger organic debris. Wipe the catheter with a wipe, focusing on the distal tip.
- Cleaning / pre-disinfection: performed in a suitable room, immerse and brush the probe in a bath of solution dedicated to the cleaning of endoscopic tools (for example: Umonium Instrument from Huckert's International).
- Disinfection: performed in a suitable room, immerse and brush the probe in a bath of solution dedicated to the cleaning of endoscopic tools (for example: Umonium Instrument from Huckert's International).

A schematic procedure for this reprocessing is presented in Appendix A of this manual.

#### 4.4.2 Automatic Reprocessing

#### Warnings and precautions for automatic reprocessing

- Ensure that the washer-disinfector is specifically designed for cleaning and disinfecting flexible endoscopic instruments (temperature sensitive above 60°C) by a manufacturer of washers-disinfectors.
- Ensure that the washer-disinfector complies with international standards ISO 15883-1 and 15883-4.
- Perform regular maintenance of the washer-disinfector, in accordance with the manufacturer's recommendations.
- Consult and follow the washer-disinfector operating instructions.
- Make sure that the *Iriscope* probe is properly supported in the washer-disinfector, for example by being placed in a basket/container provided for this purpose.
- All shut-off valves must be open.
- Do not overload the washer-disinfector.
- Remove the *Iriscope* probe from the washer-disinfector as soon as the procedure is finished as this helps to prevent the onset of corrosion.

The *Iriscope* probe can be reprocessed in any washer-disinfector, provided that it is suitable for heat-sensitive endoscopic tools (temperature below 60°C), and that it complies with international standards ISO 15883-1 and 15883-4.

Before reprocessing in the washer-disinfector, wipe the *Iriscope* probe with a damp (at least drinking water), soft, lint-free, single-use cloth. The *Iriscope* probe connector should not be wiped with this cloth, but if it is wiped, make sure that it has its cap on to protect the connections.

Then place the probe in the washer-disinfector so that it is securely held, e.g. by being placed in a basket/container provided. Ensure that the probe connector has its cap on. Start the disinfection cycle according to the manufacturer's recommendations. The following active ingredients are compatible and validated for the probe:

- Glutaraldehyde, we recommend the Olympus range: an ETD washer-disinfector, with Olympus Cleaner and Olympus Disinfectant solutions. See Appendix B for cycle details.

- Peracetic acid, we recommend the Olympus range: an ETD type washer-disinfector, with EndoDet, EndoDis and EndoAct solutions. See Appendix C for cycle details.

Drying of the probe can be carried out in the washer-disinfector if it does not exceed 60°C. If necessary, the probe can then be dried using compressed air.

If automatic drying is not possible, dry the probe with soft, lint-free, single-use cloths or suitable sterile sponges, and with compressed air if necessary.

#### 4.4.3 Hydrogen peroxide sterilisation

##### Warnings and precautions for sterilisation

- Consult and follow the instructions for use of sterilization.
- Ensure that the *Iriscope* probe is securely held in the sterilizer, by being placed in a basket/container provided for this purpose.
- Do not overload the sterilizer.
- Ensure that the correct cycle and product is used: *Sterrad 100 NX – Flex* cycle.
- Do not perform more than **20 sterilisation cycles**, the integrity of the device is not guaranteed beyond 20 cycles.
- Do not sterilise the probe in an autoclave or any other high temperature steriliser as this will destroy the probe.
- The probe must be cleaned and disinfected before being sterilised, otherwise the effectiveness of the sterilisation is not guaranteed. Refer to the methods described above (see sections 4.4.1 and 4.4.2).
- Ensure that the probe is properly dried before starting a sterilisation cycle, otherwise the effectiveness of the sterilisation is not guaranteed.

Proceed with sterilisation as follows, while respecting the manufacturer's recommendations:

1. Place the probe in a sterilisation basket (compatible with *Sterrad*).
2. Pack the basket in the dedicated sterilisation packaging and close it with the provided sticker, according to the manufacturer's (*Sterrad*) recommendations.
3. Place the probe in the *Sterrad 100 NX*, with a maximum of 1 probe per stage.
4. Start the *Sterrad 100 NX Flex* Cycle.
5. Store the probe properly, aseptically.

#### 4.5 Cleaning the connector

The electrical contacts and the ferrule must be cleaned when necessary.

##### 4.5.1 Cleaning electrical contacts

To clean connector electrical contacts, use a cotton-tip applicator or any other means that allows access to such areas to remove visible dirt.

##### 4.5.2 Cleaning the ferrule (white cylinder)

Use the tool provided with the control box to clean the ferrule, see Figure 2. The ferrule is present in the probe connector, a white cylinder, which allows good transmission of light between the control box and the probe once connected.

Figure 2 - Ferrule cleaner

Below sets out the procedure for cleaning the ferrule, illustrated in Figure 3:

1. Remove the two ferrule cleaner caps, finding a small cylindrical tip.
2. Press the tool into the cylindrical space of the connector until reaching the stop.
3. Push the tool in deep until it stops to perform cleaning. With this movement, the tip of this tool can be used to clean the ferrule tip. Perform this stage twice.
4. Now that cleaning has finished, properly replace the cap on the tool.


Figure 3 - Cleaning the ferrule

If there is insufficient light output, perform the operation again. If this situation persists, clean the control box ferrule as it may be dirty. Please see the *Iriscopes* control box manual for cleaning the ferrule.

#### 4.6 Storing the probe

After opening the packaging and between each use, the probe must be kept in a clean, well ventilated room, protected from moisture and light, and at standard temperature and pressure (see section 3.2).

The probe can be stored by being loosely rolled (curvature radius approximately 10 cm) and placed in a clean tray. If it is possible to store the probes lengthwise thus bending them less, then this is preferable.

#### 4.7 Troubleshooting

Faults	Possible causes	Solutions
No image is visible on the screen.	No probe is connected	Connect the <i>Iriscopes</i> probe.
	Problem with the control box side.	Check for possible faults in the control box user manual.
The image is blurred.	Dirty probe.	Check that the tip of the probe used is nice and clean and clean it if necessary by referring to the probe manual.
L'image est floue.	The camera shows too many signs of wear.	Check that the tip of the probe used is nice and clean and that it is not too scratched, as excessive damage to the front of the camera alters the quality of the image. Change the probe if the quality is too low for the procedure to be performed.
Problem with probe brightness.	The light emitted by the probe is of low intensity.	Detach the connector, clean the inside of the connectors – probe and control box – following the method set out, re-attach the probe to the control box. The optical fibre ferrule may be dirty.
Horizontal lines appear in the image from time to time and/or the image appears to freeze.	Signal transmission problem.	Detach the probe and check that the connector is nice and clean on both the control box side and on the probe side. Clean them if required. Re-attach the probe and check that the problem is solved. If this is not the case, the probe is damaged.

### 5 Disposal of the probe

The probe is a reusable but has a limited lifetime. The probe must be inspected before each use to ensure its good condition (see section 4.1). In case of doubt, contact Lys Medical.

If the probe is damaged, it must be disposed of in accordance with the guidelines of the institution at which it is used, and in accordance with the regulations in force.

### 6 Technical specifications

#### 6.1 Specifications

<i>Iriscopes</i> probe	<i>Iriscopes</i> 17	<i>Iriscopes</i> 18	<i>Iriscopes</i> 20
Distal tip diameter [mm]	1.7	1.85	2.0
Catheter diameter [mm]	1.5	1.5	1.5
Maximum diameter of the inserted part [mm]	1.71	1.85	2.00

<b>Iriscope probe</b>	<b>Iriscope17</b>	<b>Iriscope18</b>	<b>Iriscope20</b>
Usable length [mm]	2000	2000	2000
Minimum curvature radius [mm]	100	100	100
Field of vision [°]	120	120	120
Depth of field of vision [mm]	5 to 50	5 to 50	5 to 50
Radiance (according to EN 62471, field of vision of 100 mrad) [ $\frac{W}{m^2 sr}$ ]	109	109	109
<b>Compatible accessory</b>	<b>Iriscope17</b>	<b>Iriscope18</b>	<b>Iriscope20</b>
Minimum diameter of operating channel [mm]	1.9	2.0	2.1
<b>Storage and transport</b>	<b>Iriscope</b>		
Storage temperature [°C]	10 – 40		
Transport temperature [°C]	10 – 60		
Storage relative humidity [%]	30 – 80		
Transport relative humidity [%]	10 – 90		
Atmospheric pressure [kPa]	80 - 110		



## 6.2 Applicable standards

The *Iriscope* probe operates in accordance with European regulation *MDR 2017/745* with the following standards:

- ✓ IEC 60601-1 Medical electrical equipment – Part 1: General requirements for basic safety.
- ✓ IEC 60601-1-2: Medical electrical equipment – Part 1-2: General requirements for basic safety – Collateral standard: Electromagnetic disturbances - Requirements and tests.
- ✓ IEC 60601-2-18: Medical electrical equipment – Part 2-18: Particular requirements for the basic safety and essential performance of endoscopic equipment.
- ✓ IEC 60601-2-57: Medical electrical equipment – Part 2-57: Particular requirements for the basic safety and essential performance of non-laser light source equipment intended for therapeutic, diagnostic, monitoring and cosmetic/aesthetic use
- ✓ IEC 62417: Photobiological safety of lamps and lamp systems.
- ✓ ISO 10993-1: Biological evaluation of medical devices — Part 1: Evaluation and testing within a risk management process
- ✓ ISO 8600-1: Optics and photonic – Medical endoscopes and endotherapy devices – Part 1: General requirements
- ✓ ISO 8600-7: Optics and photonic – Medical endoscopes and endotherapy devices – Part 7: Basic requirements for medical endoscopes of water-resistant type
- ✓ ISO 17664: Processing of health care products — Information to be provided by the medical device manufacturer for the processing of medical devices

## 7 Symbols used

<b>Symbols</b>	<b>Meaning</b>
	Medical device.
	Unique device identifier.
	Catalogue reference.
	CE Marking. The device complies with European regulation <i>MDR 2017/745</i> .
	Medical device manufacturer.
	Medical device date of manufacture, in Belgium.
	Batch number.

	<p>Keep away from light.</p>
	<p>The operating instructions must be read prior to using the device.</p>
	<p>Non-sterile.</p>
	<p>Do not use if the packaging is damaged.</p>
	<p>Protect from moisture.</p>
	<p>Fragile, handle with care.</p>
	<p>Type BF applied part electrical safety.</p>
	<p>Moisture threshold: relative humidity between 30 and 80% in the operating environment.</p>
	<p>Atmospheric pressure threshold: between 800 and 1100 hPa in the operating environment.</p>
	<p>Temperature threshold: between 10 and 40°C in the operating environment.</p>

EN

## Procedure for the cleaning and disinfection of the *Iriscope* probe


### 1 Pre-cleaning

1 Wipe with a wipe soaked in quaternary ammoniums

2 Proceed from the connector to the head

3 Concentrate on joints and the camera

4 Leave to dry for 5 mins

### 2 Cleaning/Pre-disinfection

25mL/5L water

1 Use 10% quaternary ammoniums\*

2 Immerse and scrub the probe on joints

3 Rinse with potable water

4 Use 0.5% quaternary ammoniums\*

5 Immerse and scrub the probe on joints

6 Rinse with potable water

### 3 Desinfectie

100mL/1L water

1 Use 10% quaternary ammoniums\*

2 Immerse and scrub the probe on joints

3 Immerse the probe for 1h

4 Rinse with potable water

5 Dry fully and store aseptically\*\*

\* We recommend the Umonium® range by Hucker's International (Umonium® Neutralis Tissues & Umonium® Instruments)  
 \*\* If the probe is stored for more than 7 days, further disinfection is required before use

*Iriscope*

Lys Medical S.A.  
 Rue Auguste Piccard, 48  
 B-6041 Gosselies

B. Instructions for automatic reprocessing -

	Temperature	Duration	Agent	Dosage	Type of water
1) Prewashing	20°C	5 minutes	/	/	At least potable
2) Washing	48°C	5 minutes	pH neutral enzymatic product	0,6% (6 mL/L)	
3) Rinsing	20°C	5 minutes	/	/	
4) Disinfection	58°C	5 minutes	20% glutaraldehyde type product	1,2% (12 mL/L)	
5) Rinsing	20°C	5 minutes	/	/	
6) Drying	57°C	5 minutes	/	/	

We recommend *Olympus Cleaner* for washing and *Olympus Disinfectant* for disinfection.

C. Instructions for automatic reprocessing

	Temperature	Duration	Agent	Dosage	Type of water
1) Prewashing	20°C	4 minutes	/	/	At least potable
2) Washing	35°C	3 minutes	"EndoDet" Surfactant-based detergent	0,6% (6 mL/L)	
3) Disinfection	35°C	5 minutes	"EndoAct" Activator additive (phosphates + sodium hydroxide)	1,2% (12 mL/L)	
			"EndoDis" Peracetic acid	1,2% (12 mL/L)	
4) Prerinsing	20°C	4 minutes	/	/	
5) Final rinsing	20°C	4 minutes	/	/	
6) Drying	57°C	4 minutes	/	/	

We recommend Olympus products, *EndoDet* for washing and *EndoAct* + *EndoDis* for disinfection.

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