Iriscope[®]

Iriscope Guide Sheath User Manual

***Lys**
medical

Lys Medical SA
Rue Clément Ader, 10
6041 Gosseles
Belgium

info@lysmedical.com | www.lysmedical.com | +32 2 850 05 93

Version 5 – 2026-02

9. Do not use the device if the expiry date has passed. After this date, the cleanliness of the product is no longer guaranteed.

2 Description

Iriscope Guide Sheath is a catheter supporting *Iriscope* probes or any other endoscopic tool.


Figure 1: *Iriscope Guide Sheath* catheter, with (1) the access Luer connector, (2) the catheter, and (3) the optional T-shaped flush adapter.

This catheter is intended for exploring the airways for a maximum duration of 60 minutes.

Two small stopper clips are supplied with the guide. These can be used as markers to maintain the correct relative position of *Iriscope* probes and biopsy tools with the *Iriscope Guide Sheath* catheter.


Figure 2: Stopper clips for *Iriscope* probes and biopsy tools.

3 Precautions before use

3.1 Storage

⚠ Warnings and precautions

1. Be sure to comply with the storage conditions set out in this manual as improper storage could lead to an infection control risk and/or damage to the device.
2. Take care not to leave the catheter exposed to direct sunlight (UV) as this could speed up its deterioration and greatly reduce the lifetime of the catheter, or even render it unusable.

The catheter in its original packaging must be kept in a dry place away from moisture and light in a clean room, at standard temperature and pressure.

1 Important Information

These operating instructions must be read carefully prior to using the *Iriscope Guide Sheath* catheter. They contain the information required to use it safely.

This user manual describes the *Iriscope Guide Sheath* catheter, its operation and its safety precautions. It is not intended to describe clinical procedures. It is thus imperative that users have been trained in endoscopic techniques to use the *Iriscope Guide Sheath*, and be aware of all precautions, warnings, indications et contra-indications set forth in this user manual.

All serious incidents related to the *Iriscope Guide Sheath* catheter must be reported immediately to Lys Medical (vigilance@lysmedical.com) and to the responsible authority in the country where you are based.

Failure to follow the operating instructions could lead to injury to the patient and/or damage to the device. Lys Medical shall not be held liable in the event of improper use.

1.1 Intended use

The *Iriscope Guide Sheath* catheter is intended for use through bronchoscopes to extend the operating channel and thus serve as a support for compatible endoscopic tools (e.g. *Iriscope* probe, biopsy forceps, etc.). It is compatible with bronchoscopes which comply with the specifications given in section 7.1.

1.2 User qualifications

This device is intended for use by a physician specializing in bronchial endoscopy.

These operating instructions only deal with the operation of the *Iriscope Guide Sheath* catheter and its use, and endoscopic treatment and techniques are not the subject of this user manual.

1.3 Indications for use

Iriscope Guide Sheath is a single-use catheter. It must be discarded after use.

Endoscopic tools (brush, biopsy forceps, probes) can be inserted into the catheter. Tools that can be fitted must comply with the specifications given in section 7.1. However, it is recommended to test the insertion of tools in the guide beforehand (see section 3.2).

The catheter is reserved for use in adults, for the exploration of the respiratory tract, as far as the sub-segmental bronchi.

1.4 Precautions for use

It is important to heed these warnings and precautions to ensure patient safety.

Also take care to carefully read the precautions regarding catheter storage in section 3.1.

⚠ Warnings and precautions

1. The catheter must only be used by physicians trained in clinical endoscopic procedures and techniques.
2. Before use, check that the packaging is not damaged. Never use the catheter if the packaging is damaged.
3. Always check the compatibility of the *Iriscope Guide Sheath* catheter with the endoscope operating channel in advance, as well as the compatibility of endoscopic tools with *Iriscope Guide Sheath*.
4. Always preload the catheter with an endoscopic tool before use. Inserting the catheter empty could damage the catheter, endoscope, or the patient's anatomical pathways.
5. Never force the insertion or removal of the *Iriscope Guide Sheath* catheter in the endoscope operating channel.
6. Never force the insertion or removal of endoscopic tools in the *Iriscope Guide Sheath* catheter.
7. In case of malfunction during the procedure, halt the procedure and remove the *Iriscope Guide Sheath* catheter from the patient.
8. Be careful not to bend the catheter as this could permanently damage it.

3.2 Compatibility with endoscopes/endoscopic tools

It is important to check the compatibility between the *Iriscope Guide Sheath* catheter and the operating channel of the endoscope and endoscopic tools used, according to the table 7.1.

However, it cannot be guaranteed that instruments selected on the basis of the minimum or maximum recommended diameter will be compatible with the catheter. It is therefore advisable to perform a preliminary insertion test with the tools you wish to use to ensure their mechanical compatibility.

To do so, insert the catheter into the endoscope and check that insertion takes place without forcing, and that the catheter slides in correctly without catching. The use of medical lubricant can help with catheter insertion. Likewise, insert the endoscopic tools into the catheter and check that insertion takes place without forcing, and that the tools slide correctly without catching. The use of medical lubricant can help with the insertion of tools.


Incompatibility could result in damage to the *Iriscope Guide Sheath* catheter and/or injury to the patient.

4 Use


4.1 Preparation of the catheter

Before inserting the catheter into the endoscope:


1. If a biopsy tool is to be used in the *Iriscope Guide Sheath*, insert it into the channel until the end effector is at the distal exit. Place a stopper clip on the biopsy tool to mark the position of the catheter entry.

2. Screw the T-shaped flush adapter onto the *Iriscope Guide Sheath* catheter and insert the *Iriscope* probe until its tip is flush with the distal end.

3. Tighten the clamping screw of the T-shaped flush adapter to secure it on the *Iriscope* probe and seal the inlet. Flushing with a syringe makes it possible to verify the seal.

4. The T-shaped flush adapter may be left on the *Iriscope* probe to mark the position of the catheter inlet.

4.2 Inserting the catheter in the endoscope

Insert the catheter tip preloaded with the *Iriscope* probe into the endoscope operating channel. Be careful when pushing the catheter to avoid bending it, which would damage it irreversibly. The use of medical lubricant also facilitates its insertion. If resistance is encountered, do not force it – remove the catheter to check where the resistance is coming from.

4.3 Use with endoscopic tools

Once in position or during the procedure, the tools pre-loaded in the *Iriscope Guide Sheath* catheter can be removed. Different tools can be inserted or removed during the procedure as long as they are compatible with the dimensions of the guide catheter.

In case of pronounced bending of the catheter, a risk of kinking may arise. Bending the catheter to a radius of curvature of less than 15 mm is not recommended. In case of kinking, it may not be possible to insert other tools. Likewise, in the event of kinking, it may be difficult to remove tools. It is recommended to remove the guide catheter from the endoscope to de-kink the catheter in order to remove trapped tools. Dispose of the kinked *Iriscope Guide Sheath* catheter.

Take care not to force when inserting or removing endoscopic tools if resistance is encountered, as this could irreversibly damage the equipment.

The T-shaped flush accessory (Figure 3) allows flushing around the *Iriscope* probe inserted into the catheter during use. Once in position, tighten the locking ring (without forcing, to avoid damaging the device) if it is not already tightened. Use a syringe with a tip compatible with "Luer" type connectors. Check that there is no proximal flow.


Figure 3 : Tuohy-Borst adapter with (1) the clamping ring that allows the tightness around the probe and the insertion friction to be adjusted and (2) access for the syringe.

5 Troubleshooting

Faults	Possible causes	Solutions
The passage of an endoscopic tool is blocked.	The catheter is kinked. The diameter of the tool is not compatible.	Replace the catheter. Use another tool with a compatible diameter.
The "Luer" type connector has come off.	Catheter is stuck. The connector did not resist the force exerted.	Remove the endoscope to free the catheter from the patient. Replace catheter.

6 Disposal of the catheter

This is a single-use catheter. After use, dispose of it in accordance with your institution's guidelines and all applicable regulations.

7 Technical specifications

7.1 Specifications

<i>Iriscope Guide Sheath</i> catheter	GS-15-19-P	GS-16-19-B	GS-20-27-B
Maximum outer diameter of the inserted part [mm]	1,9	1,9	2,7
Nominal inner diameter [mm]	1,5	1,6	2,0
Working length [cm]	100	100	100
Minimum radius of curvature [mm]	15	15	15
Compatibility			
Maximum diameter of endoscopic tools [mm]	1,45	1,45	1,9
Minimum diameter of endoscope operating channel [mm]	2,0	2,0	2,8
Storage and transport			
Storage and operating temperature [°C]			0 - 40
Transport temperature [°C]			-20 - 60
Relative humidity during storage and operation [%]			30 - 80
Transport relative humidity [%]			10 - 90
Atmospheric pressure for storage and operation [kPa]			80 - 110
Atmospheric pressure for transport [kPa]			70 - 110

7.2 Applicable standards

The *Iriscope Guide Sheath* catheter operates in accordance with European regulation *MDR 2017/745* and with the following standards:

- ✓ ISO 10993-1: Biological evaluation of medical devices — Part 1: Evaluation and testing within a risk management process.

8 Symbols used

Symbols	Meaning
	Medical device.
	Unique device identification.
	Catalogue reference.
	CE marking. The device complies with European regulation <i>MDR 2017/745</i> .
	Medical device manufacturer.
	Batch number.
	Use-by date.
	Non-sterile device.
	Do not reuse.
	Keep out of direct sunlight.
	See instructions for use.
	Do not use if the packaging is damaged.
	Protect from moisture.
	Fragile, handle with care.
	Stopper clip supplied with guide sheath (see section 2).
	<i>Iriscope Guide Sheath</i> inner diameter.
	<i>Iriscope Guide Sheath</i> outer diameter.
	Working length.

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