Iriscope[®]

Iriscope EagleCath

User manual

Lys Medical SA

Rue Clément Ader, 10

6041 Gosselies

Belgium

info@lysmedical.com | www.lysmedical.com | +32 2 850 05 93

Version 1 – 2026-06

1 Important information

This user manual must be read carefully before using the *Iriscope EagleCath* catheter, it contains important information for safe use.

This user manual describes the *Iriscope EagleCath* catheter, its operation and its usage precautions. It is not intended to describe clinical procedures. It is therefore essential that users have been trained in endoscopic techniques in order to use *Iriscope EagleCath*, and that they are familiar with the precautions, warnings, indications and contra-indications contained in these instructions for use.

Any serious incident relating to the *Iriscope EagleCath* catheter must be reported immediately to Lys Medical (vigilance@lysmedical.com) and to the competent authority of the country in which you are established.

Failure to follow the instructions for use may result in injury to the patient and/or damage to the device. Lys Medical accepts no liability in the event of incorrect use.

1.1 Intended use

The *Iriscope EagleCath* catheter is intended to be used through bronchoscopes to extend their working channel and thus serve as a guide for compatible endoscopic tools (e.g. *Iriscope* probe, biopsy forceps, etc.). The device is single-use and must therefore be disposed of after use. It is compatible with bronchoscopes that meet the characteristics given in section 7.1.

1.2 User qualification

This device is intended to be used by a physician specialised in bronchial endoscopy.

These instructions for use deal only with the operation of the *Iriscope EagleCath* catheter and its use, endoscopic techniques and treatments are not covered by these instructions for use.

1.3 Indication for use

Iriscope EagleCath is a single-use catheter. After use, it must be disposed of in accordance with the indications given in section 6.

Endoscopic tools (biopsy forceps, probes) can be inserted into the catheter. The tools that can be inserted into it must meet the specifications given in section 7.1. The catheter is designed exclusively to insert simultaneously an *Iriscope 12* probe and a biopsy tool (forceps or brush) of 1.0 mm in diameter. The catheter can guide a single endoscopic probe or biopsy tool with a maximum external diameter of 1.50 mm (over a length of 105 cm) on its own in its channel. However, it is recommended to test the insertion of the tools into the guide beforehand (see section 3.2).

The catheter is intended for use in adults only, for the exploration of the airways, down to the subsegmental bronchi.

1.4 Precautions for use

It is important to observe these warnings and precautions to ensure patient safety.

Also be sure to carefully read the precautions concerning storage of the catheter, which are located directly in section 3.1.

Warnings and precautions

1. The catheter may only be used by physicians trained in clinical endoscopy procedures and techniques.
2. Before use, check that the packaging is not damaged. Never use the catheter if the packaging is damaged.
3. Always check beforehand the compatibility of the *Iriscope EagleCath* catheter with the working channel of the endoscope, as well as the compatibility of the endoscopic tools with *Iriscope EagleCath*.
4. Always preload the catheter with an endoscopic tool before use (at minimum an *Iriscope 12* probe). Inserting the catheter empty could damage the catheter, the endoscope, or the patient's anatomical pathways.
5. Never force the insertion or withdrawal of the *Iriscope EagleCath* catheter into the working channel of the endoscope.
6. Never force the insertion or withdrawal of the endoscopic tools into the catheter *Iriscope EagleCath*.

7. In the event of a malfunction during the procedure, discontinue the procedure and withdraw the *Iriscope EagleCath* catheter from the patient.
8. Be careful not to bend the catheter, as this could damage it permanently.
9. Do not use the device if the expiration date has passed. Beyond this date, the cleanliness of the product is no longer guaranteed.
10. Do not insert the catheter without a method to determine its position (and insertion depth) in the bronchi, in order to avoid any risk of pneumothorax.
11. Never steer the catheter while the steerable section is still in the working channel of the bronchoscope.
12. Do not apply pressure to the catheter, particularly to the steerable section, as this could damage it permanently.

2 Description

Iriscope EagleCath is a steerable catheter that can simultaneously support an *Iriscope 12* probe and a biopsy tool (forceps or brush) with an external diameter of 1.0 mm, or, on its own, any other endoscopic tool with a maximum external diameter of 1.50 mm.


Figure 1: Iriscope EagleCath catheter, with (1) the handle comprising (1a) the access "luer", (1b) the actuation lever, (2) the catheter with (2a) the steerable section, and (3) the fluidic connector allowing the leak-tight entry of two tools as well as flushing.

The catheter is supplied in packaging that guarantees its cleanliness and is ready to use for bronchoscopic procedures. This catheter is intended to explore the airways, for a duration of maximum 60 minutes.

This steerable catheter is connected to a handle that allows the catheter to be steered bidirectionally. The actuation lever (*Figure 1*, (1b)) must be pressed to unlock the steering. When the actuation lever is released, its position is locked. As the catheter transmits rotation, it is important to perform the rotation by manipulating the catheter rather than the handle alone, in order to avoid twisting the catheter.

It is supplied with a fluidic connector (*Figure 1*, (3), details in *Figure 2*) that allows a fluid to be injected to facilitate vision in the peripheral zones of the lung. A syringe with a standard luer-type tip can be connected to the flush port (see *Figure 2*). When using both tools simultaneously, the Y-connector allows the *Iriscope 12* probe to be passed (always in port (B), see *Figure 2*) and the biopsy tool (always in port (A), see *Figure 2*). Ports (A) and (B) are made leak-tight by turning the tightening screw (as shown in *Figure 2*). Tightening the entry ports also fixes the position of the tool passing through them.

When using a single tool that does not require fluid injection, the fluidic connector can be detached from the access “luer” of the handle (Figure 1, (1a)) for easier insertion into the catheter channel.


Figure 2: Fluidic connector, with the two entry ports for tools (A) and (B), as well as the flush port (indicated by a syringe symbol).

3 Precautions before use

3.1 Storage

Warnings and precautions

1. Ensure that the storage conditions mentioned in this manual are observed, as improper storage may lead to risks of infection and/or degradation of the device.
2. Be careful not to leave the catheter exposed to sunlight (UV), which could accelerate its deterioration and greatly reduce the service life of the catheter, or even make it unusable.

The catheter in its original packaging must be kept in a dry place, protected from moisture and light, in a clean room, under standard temperature and pressure.

3.2 Compatibility with endoscopes/endoscopic tools

It is important to verify the compatibility between the *Iriscope EagleCath* catheter and the working channel of the endoscope as well as the endoscopic tools used, according to the table in section 7.1.

However, it is not guaranteed that instruments selected on the basis of the recommended minimum or maximum diameter will be compatible with the catheter. It is therefore advisable to carry out a prior insertion test with the tools you wish to use to ensure mechanical compatibility between them.

To do this, insert the catheter, preloaded with at least one compatible tool, into the endoscope and check that insertion occurs smoothly without excessive force, and that the catheter moves freely inside without being blocked. Likewise, insert the endoscopic tools into the catheter and check that insertion occurs smoothly without excessive force, and that the tools move freely inside without being blocked.


An incompatibility could result in damage to the *Iriscope EagleCath* catheter and/or injury to the patient.

4 Use


4.1 Catheter preparation

Before inserting the catheter into the endoscope:

1. If the fluidic connector is not already connected to the luer of the handle and it is needed, screw it on.

2. Insert the *Iriscope* probe until the tip is flush with the distal end.

3. Tighten the tightening screw of the port (B) connector to fix it onto the *Iriscope* probe and make the entry leak-tight (Figure 2). Optionally, insert a forceps in port (A) until its tip appears in the image of the *Iriscope* probe. Tighten the screw of port (A) (even if no tool is inserted in it) to make the entry leak-tight. Flushing with a syringe allows the seal to be tested.
4. Check the correct steering of the catheter tip, in both directions. Return to the neutral position before insertion into the bronchoscope.

4.2 Inserting the catheter into the endoscope

Insert the tip of the catheter, preloaded with the *Iriscope* probe (and optionally a forceps), into the working channel of the endoscope. Be careful when pushing the catheter to avoid bending the catheter, which would damage it irreversibly. If resistance appears, do not force, and withdraw the catheter to check where the resistance is coming from.

4.3 Interaction with endoscopic tools

Once in position or during the procedure, it is possible to withdraw the tools preloaded in the *Iriscope EagleCath* catheter. Various tools can be inserted or withdrawn during the procedure, compatible with the dimensions of the guide catheter.

In the event of a pronounced bend in the catheter, a risk of kinking may appear. It is not advisable to bend the catheter to a bending radius of less than 20 mm. In the event of a kink, the insertion of new tools may be impossible. Likewise, in the event of a kink, the withdrawal of tools may be very difficult. It is recommended to withdraw the guide catheter out of the endoscope to straighten the catheter in order to remove the stuck tools. Dispose of the kinked *Iriscope EagleCath* catheter.

Be careful not to force during the withdrawal or insertion of the endoscopic tools if resistance appears, as this could irreversibly damage the equipment.

The fluidic connector (Figure 2) allows flushing around the *Iriscope* probe inserted in the catheter during use. Once in position, tighten the tightening rings of ports (A) and (B) (smoothly without excessive force, to avoid damaging the

equipment) if they are not already tightened. Use a syringe with a tip compatible with luer-type connectors. Check that the flow does not occur at the proximal level.

5 Malfunction

Malfunctions	Possible causes	Solutions
The passage of an endoscopic tool is blocked.	The catheter is kinked.	Replace the catheter.
	The diameter of the tool is not compatible.	Use another tool, with a compatible diameter.
The "luer" connector has become detached.	The catheter is stuck.	Withdraw the endoscope to release the catheter out of the patient.
	The connector did not withstand the force applied.	Replace the catheter.
The catheter cannot be steered.	The steering mechanism is damaged in the catheter.	Replace the catheter.
	The steering mechanism is damaged in the handle.	

6 Catheter disposal

The catheter is single-use, as is the fluidic connector. After use, dispose of them in accordance with the guidelines of your institution and the regulations in force.

7 Technical characteristics

7.1 Specifications

Catheter <i>Iriscope EagleCath</i>	SGS-DT-15-27
Maximum external diameter of the inserted part [mm]	2.70
Minimum internal diameter [mm]	1.60
Working length of the catheter [cm]	85
Working length of the channel [cm]	105
Minimum bending radius [mm]	20
Compatibility	
Maximum diameter of endoscopic tools inserted on their own [mm]	1.50
Maximum total of endoscopic tools inserted simultaneously [mm]	2.20
Minimum diameter of the working channel of the endoscope [mm]	2.80
Storage and transport	
Storage and operating temperature [°C]	0 - 40
Transport temperature [°C]	-20 - 60
Storage and operating relative humidity [%]	30 - 80
Transport relative humidity [%]	10 - 90
Storage and operating atmospheric pressure [kPa]	80 - 110
Transport atmospheric pressure [kPa]	70 - 110

7.2 Applied standards

The operation of the *Iriscope EagleCath* catheter complies with the European regulation *MDR 2017/745* and with the following standards:

- ✓ ISO 10993-1: Biological evaluation of medical devices – Part 1: evaluation and testing within a risk management process.

8 Symbols used

Symbols	Indication
	Medical device.
	Unique device identifier.
	Catalogue reference.
	CE marking. The device complies with the European regulation <i>MDR 2017/745</i> .
	Manufacturer of the medical device.
	Batch number.
	Use-by date.
	Non-sterile device.
	Do not reuse.
	Keep away from sunlight.
	Consult the instructions for use.
	Do not use if the packaging is damaged.
	Protect from moisture.
	Fragile, handle with care.
	Maximum diameter of tools inserted in <i>Iriscope EagleCath</i>
	External diameter of <i>Iriscope EagleCath</i>
	Working length

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